

**DR. MERCOLA®**  
HEALTHY PETS

# Tick Stick

tick removal tool

- Removes tick completely
- Quick and easy to use
- Safe for people too!



Patents  
USP 5.876.409 - EP 0 821 571  
and other patents

**DISTRIBUTED BY:** NHP  
125 SW 3rd Place  
Cape Coral, FL 33991 USA  
(877) 985-2696



1961-v111N

## Instructions for use:

1. Choose the most appropriate hook, according to the size of the tick.
2. Engage the hook by approaching the tick from the side until it is held.
3. Lift the hook very lightly and twist it until the tick detaches by itself.

